

SELF-REPORTED IMPROVEMENTS IN ENTERAL NUTRITION TOLERANCE USING MOBILITY+ ELASTOMERIC ENTERAL FEEDING SYSTEM

Edel Keaveney, PhD (Nutrition); Rachael Connolly, BSc (Human Nutrition); Warner Figueras, BEng (Biomedical Engineering). Rockfield Medical Devices Inc., c/o mHUB Chicago, 1623 W Fulton Street, Chicago, IL 60612, USA

@ info@rockfieldmd.com

www.rockfieldmd.com



BACKGROUND

- Enteral Nutrition (EN) intolerance is prevalent across care settings and may hamper EN delivery and negatively affect patients and caregivers¹, yet few studies have explored how innovative delivery systems may improve tolerance.
- This study evaluates self-reported improvements in EN feeding tolerance among early users of Mobility+ Elastomeric Enteral Feeding System.

METHODS


- Patients (n=17, 13M, 4F) were early users of the newly marketed Mobility+ (elastomeric, non-electric, portable enteral feeding system) for some, or all, their home EN feeding needs.
- Following use of Mobility+ over a minimum of two weeks, with the same EN formula previously used, patients or caregivers self-reported qualitative feedback via email and/or phone/video call.



MOBILITY+ COMPONENTS

- 1 Feeding pouch; reservoir fillable with 500ml (16.9oz) of feed
- 2 Filling-set tubing: used to fill pouch
- 3 Giving-set ('feeding-set') tubing: delivers feed to feeding tube

Patients Conditions (Age range 7-76 years)




CONCLUSIONS

- Use of Mobility+, which can be considered an alternative to simple bolus/gravity systems and more complex pump systems², was associated with self-reported improvements in EN feeding tolerance.
- This elastomeric enteral feeding system may provide a better-tolerated alternative for EN patients. Further research is warranted to explore underlying mechanisms and impact on patient outcomes.


RESULTS

- Of 17 EN patients, 13 had prior intolerance, 3 had no prior intolerance, and in one patient it was unclear.
- 13 (77%) patients self-reported improvements in EN feeding tolerance, following use of Mobility+, compared to previous systems.
- Of the 13 patients who had prior intolerance, all 13 (100%) self-reported improvements in EN feeding tolerance.

% Patients that Reported Improved Tolerance



Reported Symptom Improvements



"With my feeding pump, I felt nauseous. I was throwing up and in a lot of pain. But I haven't had any symptoms since starting Mobility+"

- Female Patient, age 41, Ohio

- Also, 3 patients reported **general tolerance improvements**.
- 4 patients reported **tolerating a higher EN feed rate** than with previous delivery systems.

REFERENCES

1. Boullata J, Carney LN, Guenter P, eds. A.S.P.E.N. Enteral Nutrition Handbook. American Society for Parenteral and Enteral Nutrition; 2010.
2. Flow Rates of a Novel Enteral Feeding System. Connolly R, Mayne D, Keaveney, E. J Parenter Enteral Nutr. 2024 FEB;48(S1): S73-S225 P27 <https://doi.org/10.1002/jpen.2604>

rockfieldmedicaldevices

rockfieldmedicaldevices

Mobility+ Enteral Feeding System

Rockfield Medical Devices