



Mobility+ Enteral Feeding System Flow Rate Guide

The Flow Rate Guide displays Mobility+ flow rates in 11 groups (A-K), from slowest to fastest. These groupings are for guidance purposes only. Actual flow rates may fall outside of the ranges indicated for each group due to a variety of factors.

FLOW	SLOW  FAST										
GROUP	A	B	C	D	E	F	G	H	I	J	K
Approx flow rate (ml/hr)	<40	41-70	71-100	101-130	131-180	181-250	251-350	351-600	601-1000	1001-1500	>1500
Approx Delivery time (hh:mm)	>12:00	12:00 - 07:00	07:00 - 05:00	05:00 - 03:45	03:45 - 02:45	02:45 - 02:00	02:00 - 01:30	01:30 - 00:45	00:45 - 00:30	00:30-00:20	<00:20

* Not suitable for use

NESTLÉ®

GIVING SET CLAMP COLOR 	WHITE	GREEN	BLUE	RED
DIABETISOURCE® AC	C	G	H	I
FIBERSOURCE® HN	D	G	H	I
ISOSOURCE® 1.5	C	E	G	H
ISOSOURCE® HN	F	H	I	J
NUTREN® 1.0	H	I	J	K
NUTREN® 1.0 FIBER	F	H	I	J
NUTREN® 1.5	C	F	G	H
NUTREN® 2.0	*	B	C	E
NUTREN® JUNIOR (VANILLA)	H	I	J	K
NUTREN® JUNIOR FIBER (VANILLA)	G	I	I	K
PEPTAMEN® 1.0	D	F	H	H
PEPTAMEN® 1.0 JR	E	G	H	I
PEPTAMEN® 1.5	B	D	E	G
PEPTAMEN® 1.5 (VANILLA)	B	D	F	G
PEPTAMEN® 1.5 JR	B	D	F	G

FLOW	SLOW FAST										
GROUP	A	B	C	D	E	F	G	H	I	J	K
Approx flow rate (ml/hr)	<40	41-70	71-100	101-130	131-180	181-250	251-350	351-600	601-1000	1001-1500	>1500
Approx Delivery time (hh:mm)	>12:00	12:00 - 07:00	07:00 - 05:00	05:00 - 03:45	03:45 - 02:45	02:45 - 02:00	02:00 - 01:30	01:30 - 00:45	00:45 - 00:30	00:30-00:20	<00:20

* Not suitable for use

NESTLÉ[®]

GIVING SET CLAMP COLOR	WHITE	GREEN	BLUE	RED
PEPTAMEN [®] INTENSE VHP	E	G	H	I
REPLETE [®]	D	G	H	I
REPLETE [®] FIBER	E	G	H	I
COMPLEAT [®] STANDARD 1.4	*	C	E	G
COMPLEAT [®] PEDIATRIC PEPTIDE 1.5	*	*	B	D
COMPLEAT [®] PEPTIDE 1.5	*	*	B	E

ABBOTT[®]

GIVING SET CLAMP COLOR	WHITE	GREEN	BLUE	RED
GLUCERNA [®] 1.5	B	C	E	F
JEVITY [®] 1.0	A	C	E	F
JEVITY [®] 1.2	C	E	F	H
JEVITY [®] 1.5	A	B	C	E
PEDIASURE [®] PEPTIDE 1.0	H	I	J	K
PEDIASURE [®] PEPTIDE 1.5	E	G	H	H
OSMOLITE [®] 1.2	E	G	H	I
OSMOLITE [®] 1.5	C	E	F	G
TWOCAL [®] HN	C	F	G	H
VITAL [®] 1.5	D	G	H	I

FLOW	SLOW FAST										
GROUP	A	B	C	D	E	F	G	H	I	J	K
Approx flow rate (ml/hr)	<40	41-70	71-100	101-130	131-180	181-250	251-350	351-600	601-1000	1001-1500	>1500
Approx Delivery time (hh:mm)	>12:00	12:00 - 07:00	07:00 - 05:00	05:00 - 03:45	03:45 - 02:45	02:45 - 02:00	02:00 - 01:30	01:30 - 00:45	00:45 - 00:30	00:30-00:20	<00:20

* Not suitable for use

KATE FARMS[®]

GIVING SET CLAMP COLOR	WHITE	GREEN	BLUE	RED
PEDIATRIC PEPTIDE 1.5	E	G	H	I
PEDIATRIC PEPTIDE 1.5 (VANILLA)	D	F	H	H
PEDIATRIC STANDARD 1.2	E	H	I	I
PEPTIDE 1.5	B	E	G	H
PEPTIDE 1.5 (VANILLA)	E	G	H	I
STANDARD 1.0	F	H	H	J
STANDARD 1.4	D	F	G	H
STANDARD 1.4 (VANILLA)	D	G	H	I

NUTRICIA[®]


GIVING SET CLAMP COLOR	WHITE	GREEN	BLUE	RED
NEOCATE [®] SPLASH	I	J	K	K

FACTORS AFFECTING FLOW RATE

For the Flow Rate Guide shown in previous pages, the Mobility+ was tested under the following controlled conditions:

- Using a single batch of each formula listed
- The formula and environmental temperature were maintained at 68-75°F (room temperature)
- The Mobility+ was positioned at stoma height
- The Mobility+ was filled with 500g of formula
- Tests were initiated immediately after filling

In everyday use conditions, several factors, listed below, may affect Mobility+ flow rates.





FORMULA VARIABILITY

Flow rates may be impacted by formula manufacturing variables with ingredients as well as formula temperature while feeding (as described below).

MOBILITY+ POUCH POSITIONING (HEIGHT)

If Mobility+ is positioned below stoma height, (e.g. worn in a fanny pack on the waist) flow rates may be slower. If Mobility+ is positioned above stoma height (e.g. in a backpack) flow rates may be faster. Mobility+ should not be worn greater than ± 6 inches above or below the stoma site.




TEMPERATURE

Environmental temperatures higher than 75°F (room temperature) may result in faster flow rates, while lower temperatures may result in slower flow rates. Mobility+ should not be used in environmental temperatures above 95°F or below 59°F. Using refrigerated formula will result in slower flow rates.

EXTERNAL PRESSURE

External pressure such as squeezing or laying on Mobility+ may result in faster flow rates.



MULTIPLE USE

Mobility+ can be re-used multiple times within 24 hours once components are flushed between use (as per instructions for use). Flow rates may be slower with multiple use.

FILL VOLUME

Do not fill Mobility+ with less than its minimum fill volume of 250ml (8.5 oz) or more than its maximum fill volume of 500ml (16.9 oz). The lower the fill volume of the Mobility+ the quicker the delivery time will be.

