

Safety & Efficacy in Tube-fed patients



A study that looked at how well Mobility+, a new elastomeric enteral feeding system, worked for adults who normally use gravity, a pump, or a syringe to get their nutrition through a feeding tube at home.

Who took part?

- ▶ 17 adults (12 men, 5 women), aged 36 to 78 years
- ▶ 14 had a G-tube, 3 had a J-tube
- ▶ They had a variety of health conditions

What did participants do?

- ▶ Participants switched to Mobility+ for 2 feedings a day for 14 days
- ▶ 16 out of 17 participants used Mobility+ for at least 2 feedings a day for 14 days
- ▶ On average, each person got 1050 ml of formula and 1657 kcals per day using Mobility+

Key findings

Safe and effective

All participants received their nutrition without any adverse events or complications.

Daily activities became easier while feeding



Moving around
the house



Walking or using
public transport



Socializing and
sleeping



Participants said Mobility+ was:



Easier to carry



Quieter



More discreet than their
previous feeding systems

Overall satisfaction and quality of life ratings stayed about the same. Some users familiar with bolus feeding noted that Mobility+ takes more effort to fill, which affected their satisfaction. 11 out of 17 participants said they would like to use Mobility+ again.

Summary

Mobility+ is a **safe and effective alternative** to traditional feeding methods. It helped participants **move more freely** and **feed more discreetly**. The system had a **positive impact on everyday life** for most users.



Visit our website
to learn more
about Mobility+